




ZZZ3 Monoclonal Antibody

Catalog No	YP-mAb-02209
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZZZ3
Immunogen	The antiserum was produced against synthesized peptide derived from human ZZZ3. AA range:341-390
Specificity	ZZZ3 Monoclonal Antibody detects endogenous levels of ZZZ3 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ZZZ3;ZZ-type zinc finger-containing protein 3
Observed Band	102kD
Calculated Molecular Weight	102kD
Cell Pathway	Nucleus .
Tissue Specificity	Adipose tissue,Brain,Colon endothelium,Lymph,Placenta,Rectum tumor,Spleen,
Function	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR., similarity:Contains 1 HTH myb-type DNA-binding domain.,similarity:Contains 1 ZZ-type zinc finger.,
Background	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR., similarity:Contains 1 HTH myb-type DNA-binding domain.,similarity:Contains 1 ZZ-type zinc finger.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using ZZZ3 Monoclonal Antibody