



KLF12 Monoclonal Antibody

Catalog No	YP-mAb-02229
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	KLF12
Immunogen	Synthesized peptide derived from the Internal region of human KLF12.
Specificity	KLF12 Monoclonal Antibody detects endogenous levels of KLF12 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	AP 2rep;AP2REP;HSPC122;KLF12;Krueppel like factor 12;Kruppel like factor 12;Krueppel-like factor 12;Transcriptional repressor AP-2rep
Observed Band	44kD
Calculated Molecular Weight	44kD
Cell Pathway	Nucleus.
Tissue Specificity	Brain,Hypothalamus,Umbilical cord blood,
Function	function:Confers strong transcriptional repression to the AP-2-alpha gene. Binds to a regulatory element (A32) in the AP-2-alpha gene promoter., similarity:Belongs to the Sp1 C2H2-type zinc-finger protein family., similarity:Contains 3 C2H2-type zinc fingers.,
Background	Activator protein-2 alpha (AP-2 alpha) is a developmentally-regulated transcription factor and important regulator of gene expression during vertebrate development and carcinogenesis. The protein encoded by this gene is a member of the Kruppel-like zinc finger protein family and can repress expression of the AP-2 alpha gene by binding to a specific site in the AP-2 alpha gene promoter. Repression by the encoded protein requires binding with a corepressor, CtBP1. Two transcript variants encoding different isoforms have



been found for this gene. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using KLF12 Monoclonal Antibody