




CGR19 Monoclonal Antibody

Catalog No	YP-mAb-02242
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CGRRF1
Immunogen	The antiserum was produced against synthesized peptide derived from the C-terminal region of human CGRRF1. AA range:283-332
Specificity	CGR19 Monoclonal Antibody detects endogenous levels of CGR19 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CGR19;CGRRF1;RING finger protein 197;RNF197;Cell growth regulator with RING finger domain protein 1;Cell growth regulatory gene 19 protein
Observed Band	38kD
Calculated Molecular Weight	38kD
Cell Pathway	Nucleus . Endoplasmic reticulum .
Tissue Specificity	Ubiquitously expressed with high expression in testis and the cerebellum.
Function	function:Able to inhibit growth in several cell lines.,similarity:Contains 1 RING-type zinc finger.,
Background	function:Able to inhibit growth in several cell lines.,similarity:Contains 1 RING-type zinc finger.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using CGR19 Monoclonal Antibody