



ABHD11 Monoclonal Antibody

Catalog No	YP-mAb-02450
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ABHD11
Immunogen	The antiserum was produced against synthesized peptide derived from human ABHD11. AA range:161-210
Specificity	ABHD11 Monoclonal Antibody detects endogenous levels of ABHD11 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	sn-1-specific diacylglycerol lipase ABHD11;PP1226;EC:3.1.1.116;Alpha/beta hydrolase domain-containing protein 11;Abhydrolase domain-containing protein 11;ABHD11;WBSCR21;Williams-Beuren syndrome chromosomal region 21 protein
Observed Band	32kD
Calculated Molecular Weight	32kD
Cell Pathway	mitochondrion,
Tissue Specificity	Ubiquitously expressed.
Function	caution:It is uncertain whether Met-1 or Met-10 is the initiator.,sequence caution:Wrong choice of frame.,similarity:Belongs to the AB hydrolase superfamily.,tissue specificity:Ubiquitously expressed.,
Background	This gene encodes a protein containing an alpha/beta hydrolase fold domain. This gene is deleted in Williams syndrome, a multisystem developmental disorder caused by the deletion of contiguous genes at 7q11.23. [provided by RefSeq, Mar 2016],



matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using ABHD11 Monoclonal Antibody