



ACOT4 Monoclonal Antibody

Catalog No	YP-mAb-02471
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ACOT4
Immunogen	The antiserum was produced against synthesized peptide derived from human ACOT4. AA range:361-410
Specificity	ACOT4 Monoclonal Antibody detects endogenous levels of ACOT4 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ACOT4;acyl CoA thioesterase 4;Acyl coenzyme A thioesterase 4;PTE 2b;PTE 1b; PTE1B;PTE2B;PTE1B;Acyl-coenzyme A thioesterase 4;Acyl-CoA thioesterase 4; PTE-2b;Peroxisomal acyl coenzyme A thioester hydrolase 1b;Peroxisomal long-chain acyl-CoA thioesterase 1b;PTE-1b
Observed Band	48kD
Calculated Molecular Weight	48kD
Cell Pathway	Peroxisome .
Tissue Specificity	Strongest expression in liver and kidney and weaker expression in placenta, heart, and muscle.
Function	catalytic activity:Palmitoyl-CoA + H(2)O = CoA + palmitate.,function:Acyl-CoA thioesterases are a group of enzymes that catalyze the hydrolysis of acyl-CoAs to the free fatty acid and coenzyme A (CoASH), providing the potential to regulate intracellular levels of acyl-CoAs, free fatty acids and CoASH., similarity:Belongs to the C/M/P thioester hydrolase family.,



Background

catalytic activity: Palmitoyl-CoA + H₂O = CoA + palmitate., function: Acyl-CoA thioesterases are a group of enzymes that catalyze the hydrolysis of acyl-CoAs to the free fatty acid and coenzyme A (CoASH), providing the potential to regulate intracellular levels of acyl-CoAs, free fatty acids and CoASH., similarity: Belongs to the C/M/P thioester hydrolase family.,


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using ACOT4 Monoclonal Antibody