



CA VA Monoclonal Antibody

Catalog No	YP-mAb-02513
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CA5A
Immunogen	The antiserum was produced against synthesized peptide derived from human CA5A. AA range:171-220
Specificity	CA VA Monoclonal Antibody detects endogenous levels of CA VA protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CAV;Carbonic anhydrase VA;Carbonic anhydrase 5A;mitochondrial;Carbonate dehydratase VA;CA5;CA VA;CA5A;Carbonic anhydrase 5A mitochondrial;CA-VA
Observed Band	35kD
Calculated Molecular Weight	35kD
Cell Pathway	Mitochondrion .
Tissue Specificity	Liver,Umbilical cord,
Function	catalytic activity:H(2)CO(3) = CO(2) + H(2)O.,cofactor:Zinc.,function:Reversible hydration of carbon dioxide.,similarity:Belongs to the alpha-carbonic anhydrase family.,
Background	Carbonic anhydrases (CAs) are a large family of zinc metalloenzymes that catalyze the reversible hydration of carbon dioxide. They participate in a variety of biological processes, including respiration, calcification, acid-base balance, bone resorption, and the formation of aqueous humor, cerebrospinal fluid, saliva, and gastric acid. They show extensive diversity in tissue distribution and in their subcellular localization. CA VA is localized in the mitochondria and expressed primarily in the liver. It may play an important role in ureagenesis and



gluconeogenesis. CA5A gene maps to chromosome 16q24.3 and an unprocessed pseudogene has been assigned to 16p12-p11.2. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using CA VA Monoclonal Antibody