



CYP24A1 Monoclonal Antibody

Catalog No	YP-mAb-02570
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CYP24A1
Immunogen	The antiserum was produced against synthesized peptide derived from human Cytochrome P450 24A1. AA range:411-460
Specificity	CYP24A1 Monoclonal Antibody detects endogenous levels of CYP24A1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Cytochrome P450 24A1;CYP24;CP24;24-OHase;24 OHase;CYP24A1;1;25-dihydroxyvitamin D(3) 24-hydroxylase mitochondrial;25-dihydroxyvitamin D(3) 24-hydroxylase;mitochondrial;Vitamin D(3) 24-hydroxylase;Cytochrome P450-CC24; mitochondrial
Observed Band	38kD
Calculated Molecular Weight	38kD
Cell Pathway	Mitochondrion .
Tissue Specificity	PCR rescued clones,
Function	catalytic activity:Calcidiol + NADPH + O(2) = secalciferol + NADP(+) + H(2)O., catalytic activity:Calcitriol + NADPH + O(2) = calcitetrol + NADP(+) + H(2)O., cofactor:Heme group.,function:Has a role in maintaining calcium homeostasis. Catalyzes the NADPH-dependent 24-hydroxylation of calcidiol (25-hydroxyvitamin D(3)) and calcitriol (1-alpha,25-dihydroxyvitamin D(3)). The enzyme can perform up to 6 rounds of hydroxylation of calcitriol leading to calcitroic acid. It also shows 23-hydroxylating activity leading to 1-alpha,25-dihydroxyvitamin D(3)-26,23-lactone as end product.,similarity:Belongs to the cytochrome P450 family.,



Background

cytochrome P450 family 24 subfamily A member 1(CYP24A1) Homo sapiens
This gene encodes a member of the cytochrome P450 superfamily of enzymes. The cytochrome P450 proteins are monooxygenases which catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids. This mitochondrial protein initiates the degradation of 1,25-dihydroxyvitamin D3, the physiologically active form of vitamin D3, by hydroxylation of the side chain. In regulating the level of vitamin D3, this enzyme plays a role in calcium homeostasis and the vitamin D endocrine system. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using CYP24A1 Monoclonal Antibody