



DHODH Monoclonal Antibody

Catalog No	YP-mAb-02616
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	DHODH
Immunogen	The antiserum was produced against synthesized peptide derived from human DHODH. AA range:211-260
Specificity	DHODH Monoclonal Antibody detects endogenous levels of DHODH protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PYRD;DHODEHase;dihydroorotate dehydrogenase;Dihydroorotate dehydrogenase (quinone);mitochondrial;Dihydroorotate oxidase;DHODH; quinone)
Observed Band	48kD
Calculated Molecular Weight	48kD
Cell Pathway	Mitochondrion inner membrane ; Single-pass membrane protein .
Tissue Specificity	Skin,Testis,
Function	catalytic activity:(S)-dihydroorotate + O(2) = orotate + H(2)O(2).,cofactor: Binds 1 FAD per subunit.,pathway:Pyrimidine metabolism; UMP biosynthesis via de novo pathway; UMP from HCO(3)(-): step 4/6.,similarity:Belongs to the dihydroorotate dehydrogenase family.,
Background	The protein encoded by this gene catalyzes the fourth enzymatic step, the ubiquinone-mediated oxidation of dihydroorotate to orotate, in de novo pyrimidine biosynthesis. This protein is a mitochondrial protein located on the outer surface of the inner mitochondrial membrane. [provided by RefSeq, Jul 2008],


**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using DHODH Monoclonal Antibody