



DHRS2 Monoclonal Antibody

Catalog No	YP-mAb-02619
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	DHRS2
Immunogen	The antiserum was produced against synthesized peptide derived from human DHRS2. AA range:111-160
Specificity	DHRS2 Monoclonal Antibody detects endogenous levels of DHRS2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DHRS2;Dicarbonyl reductase HEP27;HEP27;Protein D;SDR25C1;Dehydrogenase/reductase SDR family member 2
Observed Band	27kD
Calculated Molecular Weight	27kD
Cell Pathway	Mitochondrion matrix. Nucleus. A minor fraction of the protein is translocated from the mitochondria to the nucleus, after cleavage of the targeting signal.
Tissue Specificity	Widely expressed, with highest levels in liver and kidney, followed by heart, spleen, skeletal muscle and placenta. In hemopoietic cells, expressed in dendritic cells, but not in monocytes, macrophages, granulocytes, nor in B and T lymphocytes.
Function	function:May inhibit cell replication either by catalyzing the oxidation of estrogen and androgen or by converting cortisone in cortisol.,similarity:Belongs to the short-chain dehydrogenases/reductases (SDR) family.,
Background	function:May inhibit cell replication either by catalyzing the oxidation of estrogen and androgen or by converting cortisone in cortisol.,similarity:Belongs to the short-chain dehydrogenases/reductases (SDR) family.,




matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using DHRS2 Monoclonal Antibody