




GUSBP1 Monoclonal Antibody

Catalog No	YP-mAb-02650
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GUSBP1
Immunogen	The antiserum was produced against synthesized peptide derived from human GUSBL1. AA range:51-100
Specificity	GUSBP1 Monoclonal Antibody detects endogenous levels of GUSBP1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GUSBP1;Putative beta-glucuronidase-like protein SMA3;SMA3
Observed Band	18kD
Calculated Molecular Weight	18kD
Cell Pathway	
Tissue Specificity	Ubiquitous.
Function	skeletal system development,
Background	Binds GTP, GDP and GMP.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using GUSBP1 Monoclonal Antibody