



MMP-19 Monoclonal Antibody

Catalog No	YP-mAb-02681
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MMP19
Immunogen	The antiserum was produced against synthesized peptide derived from human MMP-19. AA range:11-60
Specificity	MMP-19 Monoclonal Antibody detects endogenous levels of MMP-19 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EC:3.4.24.-;matrix metalloproteinase 19;Matrix metalloproteinase 18;Matrix metalloproteinase 19;Matrix metalloproteinase-18;MMP-19;MMP19;Matrix metalloproteinase-19;MMP18;RASI;Matrix metalloproteinase RASI;MMP-18
Observed Band	57kD
Calculated Molecular Weight	57kD
Cell Pathway	Secreted, extracellular space, extracellular matrix .
Tissue Specificity	Expressed in mammary gland, placenta, lung, pancreas, ovary, small intestine, spleen, thymus, prostate, testis colon, heart and blood vessel walls. Not detected in brain and peripheral blood leukocytes. Also expressed in the synovial fluid of normal and rheumatoid patients (PubMed:8920941).
Function	alternative products:Additional isoforms seem to exist,catalytic activity:Cleaves aggrecan at the 360-Ser- -Phe-361 site.,cofactor:Bounds 1 zinc ion per subunit.,cofactor:Calcium.,disease:May play a role in pathological processes participating in rheumatoid arthritis (RA)-associated joint tissue destruction. Autoantigen anti-MMP19 are frequent in RA patients.,domain:The conserved cysteine present in the cysteine-switch motif binds the catalytic zinc ion, thus inhibiting the enzyme. The dissociation of the cysteine from the zinc ion upon



the activation-peptide release activates the enzyme.,enzyme regulation:Strongly inhibited by TIMP-2, TIMP-3 and TIMP-4, while TIMP-1 is less efficient., function:Endopeptidase that degrades various components of the extracellular matrix, such as aggrecan and cartilage oligomeric matrix protein (comp), during development, haemostasis and pathological conditi

Background

This gene encodes a member of a family of proteins that are involved in the breakdown of extracellular matrix in normal physiological processes, such as embryonic development, reproduction, and tissue remodeling, as well as in disease processes, such as arthritis and metastasis. The encoded protein is secreted as an inactive proprotein, which is activated upon cleavage by extracellular proteases. Alternative splicing results in multiple transcript variants for this gene. [provided by RefSeq, Jan 2013],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using MMP-19 Monoclonal Antibody