




NDUFA4L2 Monoclonal Antibody

Catalog No	YP-mAb-02699
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	NDUFA4L2
Immunogen	The antiserum was produced against synthesized peptide derived from human NDUFA4L2. AA range:38-87
Specificity	NDUFA4L2 Monoclonal Antibody detects endogenous levels of NDUFA4L2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NUOMS;NADH-ubiquinone oxidoreductase MLRQ subunit homolog;NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 4-like 2;NDUFA4L2
Observed Band	15kD
Calculated Molecular Weight	15kD
Cell Pathway	mitochondrial respiratory chain complex IV,integral component of membrane,
Tissue Specificity	Adrenal gland,Muscle,
Function	similarity:Belongs to the complex I NDUFA4 subunit family.,
Background	similarity:Belongs to the complex I NDUFA4 subunit family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using NDUFA4L2 Monoclonal Antibody