



NDUFC2 Monoclonal Antibody

Catalog No	YP-mAb-02706
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	NDUFC2
Immunogen	The antiserum was produced against synthesized peptide derived from human NDUC2. AA range:51-100
Specificity	NDUFC2 Polyclonal Antibody detects endogenous levels of NDUFC2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB1:500 - 1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HLC 1;Complex I-B14.5b;Complex I B14.5b;CI-B14.5b;CI B14.5b;NDUFC2; NADH dehydrogenase [ubiquinone] 1 subunit C2;HLC1;Human lung cancer oncogene 1 protein;HLC-1;NADH-ubiquinone oxidoreductase subunit B14.5b
Observed Band	14kDa
Calculated Molecular Weight	14kDa
Cell Pathway	Mitochondrion inner membrane ; Single-pass membrane protein ; Matrix side .
Tissue Specificity	Kidney,Ovary,Umbilical cord blood,
Function	function:Accessory subunit of the mitochondrial membrane respiratory chain NADH dehydrogenase (Complex I), that is believed to be not involved in catalysis. Complex I functions in the transfer of electrons from NADH to the respiratory chain. The immediate electron acceptor for the enzyme is believed to be ubiquinone.,similarity:Belongs to the complex I NDUFC2 subunit family., subunit:Complex I is composed of 45 different subunits.,
Background	function:Accessory subunit of the mitochondrial membrane respiratory chain NADH dehydrogenase (Complex I), that is believed to be not involved in catalysis. Complex I functions in the transfer of electrons from NADH to the

respiratory chain. The immediate electron acceptor for the enzyme is believed to be ubiquinone.,similarity:Belongs to the complex I NDUFC2 subunit family., subunit:Complex I is composed of 45 different subunits.,

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images