




NICE4 Monoclonal Antibody

Catalog No	YP-mAb-02717
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	UBAP2L
Immunogen	The antiserum was produced against synthesized peptide derived from human UBAP2L. AA range:561-610
Specificity	NICE4 Monoclonal Antibody detects endogenous levels of NICE4 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	1E3D1;KIAA0144;NICE 4;NICE4;Protein NICE 4;UBAP2L;Ubiquitin-associated protein 2-like;Protein NICE-4
Observed Band	115kD
Calculated Molecular Weight	115kD
Cell Pathway	Cytoplasm, Stress granule .
Tissue Specificity	Ubiquitous.
Function	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR., similarity:Contains 1 UBA domain.,tissue specificity:Ubiquitous.,
Background	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR., similarity:Contains 1 UBA domain.,tissue specificity:Ubiquitous.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using NICE4 Monoclonal Antibody