



PDE10A Monoclonal Antibody

Catalog No	YP-mAb-02739
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PDE10A
Immunogen	The antiserum was produced against synthesized peptide derived from human PDE10A. AA range:21-70
Specificity	PDE10A Monoclonal Antibody detects endogenous levels of PDE10A protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	phosphodiesterase 10A;HSPDE10A;EC:3.1.4.17;cAMP and cAMP-inhibited cGMP 3',5'-cyclic phosphodiesterase 10A;PDE10A
Observed Band	75kD
Calculated Molecular Weight	88kD
Cell Pathway	Cytoplasm, cytosol .
Tissue Specificity	Abundant in the putamen and caudate nucleus regions of brain and testis, moderately expressed in the thyroid gland, pituitary gland, thalamus and cerebellum.
Function	alternative products:Isoforms differ in their N-terminal region,catalytic activity:Guanosine 3',5'-cyclic phosphate + H(2)O = guanosine 5'-phosphate., catalytic activity:Nucleoside 3',5'-cyclic phosphate + H(2)O = nucleoside 5'-phosphate.,cofactor:Binds 1 magnesium ion.,cofactor:Binds 1 zinc ion., domain:Composed of a C-terminal catalytic domain containing two divalent metal sites and an N-terminal regulatory domain which contains one cyclic nucleotide-binding region.,domain:The tandem GAF domains bind cAMP, and regulate enzyme activity. The binding of cAMP stimulates enzyme activity., enzyme regulation:Inhibited by dipyrindamole and moderately by IBMX. cAMP



acts as an allosteric activator.,function:Plays a role in signal transduction by regulating the intracellular concentration of cyclic nucleotides. Can hydrolyze both cAMP and cGMP, but has higher affinity for cAMP and is more efficient w

Background

The protein encoded by this gene belongs to the cyclic nucleotide phosphodiesterase family. It plays a role in signal transduction by regulating the intracellular concentration of cyclic nucleotides. This protein can hydrolyze both cAMP and cGMP to the corresponding nucleoside 5' monophosphate, but has higher affinity for cAMP, and is more efficient with cAMP as substrate. Alternatively spliced transcript variants have been described for this gene. [provided by RefSeq, Dec 2011],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using PDE10A Monoclonal Antibody