




PFAAP5 Monoclonal Antibody

Catalog No	YP-mAb-02748
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	N4BP2L2
Immunogen	The antiserum was produced against synthesized peptide derived from human N4BP2L2. AA range:358-407
Specificity	PFAAP5 Monoclonal Antibody detects endogenous levels of PFAAP5 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	92M18.3;CG005;N4BP2L2;NEDD4 binding protein 2 like 2;PFAAP5;NEDD4-binding protein 2-like 2;Phosphonoformate immuno-associated protein 5
Observed Band	67kD
Calculated Molecular Weight	67kD
Cell Pathway	nucleus,transcriptional repressor complex,extracellular exosome,
Tissue Specificity	Colon,Thalamus,
Function	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,
Background	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using PFAAP5 Monoclonal Antibody