



SENP5 Monoclonal Antibody

Catalog No	YP-mAb-02781
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SENP5
Immunogen	The antiserum was produced against synthesized peptide derived from human SENP5. AA range:651-700
Specificity	SENP5 Monoclonal Antibody detects endogenous levels of SENP5 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SENP5;Sentrin-specific protease 5;Sentrin/SUMO-specific protease SENP5; EC:3.4.22.;FKSG45
Observed Band	86kD
Calculated Molecular Weight	86kD
Cell Pathway	Nucleus, nucleolus .
Tissue Specificity	Placenta,Testis,
Function	function:Protease that catalyzes two essential functions in the SUMO pathway: processing of full-length SUMO3 to its mature form and deconjugation of SUMO2 and SUMO3 from targeted proteins. Has weak proteolytic activity against full-length SUMO1 or SUMO1 conjugates. Required for cell division., similarity:Belongs to the peptidase C48 family.,
Background	The reversible posttranslational modification of proteins by the addition of small ubiquitin-like SUMO proteins (see SUMO1; MIM 601912) is required for numerous biologic processes. SUMO-specific proteases, such as SENP5, are responsible for the initial processing of SUMO precursors to generate a C-terminal diglycine motif required for the conjugation reaction. They also have



isopeptidase activity for the removal of SUMO from high molecular mass SUMO conjugates (Di Bacco et al., 2006 [PubMed 16738315]).[supplied by OMIM, Jun 2009],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using SENP5 Monoclonal Antibody