



SENP7 Monoclonal Antibody

Catalog No	YP-mAb-02783
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SENP7
Immunogen	The antiserum was produced against synthesized peptide derived from human SENP7. AA range:991-1040
Specificity	SENP7 Monoclonal Antibody detects endogenous levels of SENP7 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SENP7;Sentrin-specific protease 7;KIAA1707;SSP2;SUSP2;SUMO-1-specific protease 2;Sentrin/SUMO-specific protease SENP7
Observed Band	112kD
Calculated Molecular Weight	112kD
Cell Pathway	Cytoplasm .
Tissue Specificity	Brain,Epithelium,Fetal kidney,Placenta,
Function	alternative products:Experimental confirmation may be lacking for some isoforms,caution:The sequence shown here is derived from an Ensembl automatic analysis pipeline and should be considered as preliminary data., function:Protease that deconjugates SUMO2 and SUMO3 from targeted proteins, but not SUMO1. Catalyzes the deconjugation of poly-SUMO2 and poly-SUMO3 chains. Has very low efficiency in processing full-length SUMO proteins to their mature forms.,similarity:Belongs to the peptidase C48 family.,
Background	The reversible posttranslational modification of proteins by the addition of small ubiquitin-like SUMO proteins (see SUMO1; MIM 601912) is required for many cellular processes. SUMO-specific proteases, such as SENP7, process SUMO



precursors to generate a C-terminal diglycine motif required for the conjugation reaction. They also display isopeptidase activity for deconjugation of SUMO-conjugated substrates (Lima and Reverter, 2008 [PubMed 18799455]).[supplied by OMIM, Jun 2009],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using SENP7 Monoclonal Antibody