



TPH2 Monoclonal Antibody

Catalog No	YP-mAb-02803
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TPH2
Immunogen	The antiserum was produced against synthesized peptide derived from human TPH2. AA range:4-53
Specificity	TPH2 Monoclonal Antibody detects endogenous levels of TPH2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EC:1.14.16.4;NTPH;TPH 2;Tryptophan 5-hydroxylase 2;Tryptophan 5-monooxygenase 2;TPH2;Neuronal tryptophan hydroxylase
Observed Band	56kD
Calculated Molecular Weight	56kD
Cell Pathway	cytosol,neuron projection,
Tissue Specificity	Brain specific.
Function	catalytic activity:L-tryptophan + tetrahydrobiopterin + O(2) = 5-hydroxy-L-tryptophan + 4a-hydroxytetrahydrobiopterin.,cofactor:Fe(2+) ion., disease:Genetic variation in TPH2 may influence susceptibility to major depressive disorder (MDD) [MIM:608516].,pathway:Aromatic compound metabolism; serotonin biosynthesis; serotonin from L-tryptophan: step 1/2., similarity:Belongs to the biopterin-dependent aromatic amino acid hydroxylase family.,similarity:Contains 1 ACT domain.,tissue specificity:Brain specific.,
Background	This gene encodes a member of the pterin-dependent aromatic acid hydroxylase family. The encoded protein catalyzes the first and rate limiting step in the biosynthesis of serotonin, an important hormone and neurotransmitter.



Mutations in this gene may be associated with psychiatric diseases such as bipolar affective disorder and major depression. [provided by RefSeq, Feb 2016],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using TPH2 Monoclonal Antibody