



USP19 Monoclonal Antibody

Catalog No	YP-mAb-02830
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	USP19
Immunogen	The antiserum was produced against synthesized peptide derived from human USP19. AA range:391-440
Specificity	USP19 Monoclonal Antibody detects endogenous levels of USP19 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EC:3.4.19.12;KIAA0891;Ubiquitin carboxyl-terminal hydrolase 19;Ubiquitin thioesterase 19;Ubiquitin-specific-processing protease 19;USP19;ZMYND9;Deubiquitinating enzyme 19;Zinc finger MYND domain-containing protein 9
Observed Band	130kD
Calculated Molecular Weight	130kD
Cell Pathway	Endoplasmic reticulum membrane ; Single-pass membrane protein .
Tissue Specificity	Brain,Epithelium,Skin,Testis,
Function	catalytic activity:Ubiquitin C-terminal thioester + H(2)O = ubiquitin + a thiol., function:Involved in the ubiquitin-dependent proteolytic pathway in conjunction with the 26S proteasome.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,sequence caution:Translation N-terminally shortened., similarity:Belongs to the peptidase C19 family.,similarity:Contains 1 MYND-type zinc finger.,similarity:Contains 2 CS domains.,
Background	catalytic activity:Ubiquitin C-terminal thioester + H(2)O = ubiquitin + a thiol., function:Involved in the ubiquitin-dependent proteolytic pathway in conjunction with the 26S proteasome.,PTM:Phosphorylated upon DNA damage, probably by



ATM or ATR.,sequence caution:Translation N-terminally shortened.,
similarity:Belongs to the peptidase C19 family.,similarity:Contains 1 MYND-type
zinc finger.,similarity:Contains 2 CS domains.,


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using USP19 Monoclonal Antibody