



USP30 Monoclonal Antibody

Catalog No	YP-mAb-02832
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	USP30
Immunogen	The antiserum was produced against synthesized peptide derived from human USP30. AA range:31-80
Specificity	USP30 Monoclonal Antibody detects endogenous levels of USP30 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Ubiquitin-specific-processing protease 30;Ubiquitin thiolesterase 30;Ubiquitin thioesterase 30;Ubiquitin carboxyl-terminal hydrolase 30;EC:3.4.19.12;USP30; Deubiquitinating enzyme 30;Ub-specific protease 30
Observed Band	60kd unspecial bands(65kd 105kd 130kd)in Hela and HCT116 cell.
Calculated Molecular Weight	60kd,unspecial bands(65kd,105kd,130kd)in Hela and HCT116 cell.
Cell Pathway	Mitochondrion outer membrane .
Tissue Specificity	Expressed in skeletal muscle, pancreas, liver and kidney.
Function	catalytic activity:Ubiquitin C-terminal thioester + H(2)O = ubiquitin + a thiol., function:May participate in the maintenance of mitochondrial morphology., similarity:Belongs to the peptidase C19 family.,tissue specificity:Expressed in skeletal muscle, pancreas, liver and kidney.,
Background	USP30, a member of the ubiquitin-specific protease family (see USP1, MIM 603478), is a novel mitochondrial deubiquitinating (DUB) enzyme (Nakamura and Hirose, 2008 [PubMed 18287522]).[supplied by OMIM, Dec 2008],




matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using USP30 Monoclonal Antibody