




ZADH2 Monoclonal Antibody

Catalog No	YP-mAb-02844
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZADH2
Immunogen	The antiserum was produced against synthesized peptide derived from human ZADH2. AA range:191-240
Specificity	ZADH2 Monoclonal Antibody detects endogenous levels of ZADH2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ZADH2;Zinc-binding alcohol dehydrogenase domain-containing protein 2
Observed Band	45kD
Calculated Molecular Weight	45kD
Cell Pathway	Peroxisome .
Tissue Specificity	Brain,Testis,
Function	similarity:Belongs to the zinc-containing alcohol dehydrogenase family. Quinone oxidoreductase subfamily.,
Background	similarity:Belongs to the zinc-containing alcohol dehydrogenase family. Quinone oxidoreductase subfamily.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using ZADH2 Monoclonal Antibody