



Amylase Monoclonal Antibody

Catalog No	YP-mAb-02858
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	AMY1A/AMY1B/AMY1C/AMY2A/AMY2B
Immunogen	The antiserum was produced against synthesized peptide derived from the N-terminal region of human AMY1/2. AA range:61-110
Specificity	Amylase Monoclonal Antibody detects endogenous levels of Amylase protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Amylase;AMY1A/AMY1B/AMY1C/AMY2A/AMY2B;Alpha-amylase 1/Pancreatic alpha-amylase/Alpha-amylase 2B;AMY1A;AMY1;AMY1B;AMY1C;Alpha-amylase 1;1;4-alpha-D-glucan glucanohydrolase 1;Salivary alpha-amylase;AMY2A;Pancreatic alpha-amylase;PA;4-alpha-D-glucan glucanohydrolase;AMY2B;Alpha-amylase 2B;4-alpha-D-glucan glucanohydrolase 2B;Carcinoid alpha-amylase
Observed Band	60kD
Calculated Molecular Weight	60kD
Cell Pathway	extracellular space,extracellular exosome,
Tissue Specificity	Glandular pool- thyroid,Saliva,Thyroid,
Function	catalytic activity:Endohydrolysis of (1->4)-alpha-D-glucosidic linkages in oligosaccharides and polysaccharides.,cofactor:Binds 1 calcium ion per subunit.,cofactor:Binds 1 chloride ion per subunit.,online information:Amylase entry,similarity:Belongs to the glycosyl hydrolase 13 family.,subunit:Monomer.,



Background

Amylases are secreted proteins that hydrolyze 1,4-alpha-glucoside bonds in oligosaccharides and polysaccharides, and thus catalyze the first step in digestion of dietary starch and glycogen. The human genome has a cluster of several amylase genes that are expressed at high levels in either salivary gland or pancreas. This gene encodes an amylase isoenzyme produced by the salivary gland. Alternative splicing results in multiple transcript variants encoding the same protein. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using Amylase Monoclonal Antibody