




IFITM2 Monoclonal Antibody

Catalog No	YP-mAb-02878
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	IFITM2
Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human IFITM2. AA range:1-50
Specificity	IFITM2 Monoclonal Antibody detects endogenous levels of IFITM2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	1 8D;18D;DSPA2c;IFITM 2;Interferon-induced transmembrane protein 2;IFITM2; Dispanin subfamily A member 2c;Interferon-inducible protein 1-8D
Observed Band	15kD
Calculated Molecular Weight	15kD
Cell Pathway	Cell membrane ; Single-pass type II membrane protein . Lysosome membrane ; Single-pass type II membrane protein . Late endosome membrane ; Single-pass type II membrane protein .
Tissue Specificity	Lung,
Function	induction:By alpha and gamma interferons.,similarity:Belongs to the CD225 family.,
Background	induction:By alpha and gamma interferons.,similarity:Belongs to the CD225 family.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using IFITM2 Monoclonal Antibody