



ARALAR Monoclonal Antibody

Catalog No	YP-mAb-03086
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SLC25A12
Immunogen	The antiserum was produced against synthesized peptide derived from human CMC1. AA range:391-440
Specificity	ARALAR Monoclonal Antibody detects endogenous levels of ARALAR protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	AGC1;Araceli hiperlarga;Aralar;Aralar1;Electrogenic aspartate/glutamate antiporter SLC25A12;mitochondrial;SLC25A12;Calcium-binding mitochondrial carrier protein Aralar1;Mitochondrial aspartate glutamate carrier 1;Solute carrier family 25 member 12
Observed Band	75kD
Calculated Molecular Weight	75kD
Cell Pathway	Mitochondrion inner membrane ; Multi-pass membrane protein .
Tissue Specificity	Expressed predominantly in the heart and skeletal muscle, weakly in brain and kidney.
Function	function:Calcium-dependent mitochondrial aspartate and glutamate carrier. May have a function in the urea cycle.,miscellaneous:Binds calcium.,similarity:Belongs to the mitochondrial carrier family.,similarity:Contains 3 Solcar repeats.,similarity:Contains 4 EF-hand domains.,tissue specificity:High levels in heart and skeletal muscle, low in brain and very low in kidney.,



Background

This gene encodes a calcium-binding mitochondrial carrier protein. The encoded protein localizes to the mitochondria and is involved in the exchange of aspartate for glutamate across the inner mitochondrial membrane. Polymorphisms in this gene may be associated with autism, and mutations in this gene may also be a cause of global cerebral hypomyelination. Alternatively spliced transcript variants have been observed for this gene. [provided by RefSeq, Apr 2012],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using ARALAR Monoclonal Antibody