




CaBP7 Monoclonal Antibody

Catalog No	YP-mAb-03089
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CABP7
Immunogen	The antiserum was produced against synthesized peptide derived from human CABP7. AA range:1-50
Specificity	CaBP7 Monoclonal Antibody detects endogenous levels of CaBP7 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Calcium-binding protein 7;CALN2;Calneuron II;Calneuron-2;CaBP7
Observed Band	25kD
Calculated Molecular Weight	25kD
Cell Pathway	Golgi apparatus, trans-Golgi network membrane ; Single-pass type IV membrane protein . Cytoplasm, perinuclear region . Cell membrane ; Single-pass type IV membrane protein .
Tissue Specificity	Brain,Epithelium,
Function	similarity:Contains 2 EF-hand domains.,
Background	similarity:Contains 2 EF-hand domains.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using CaBP7 Monoclonal Antibody