



CIB3 Monoclonal Antibody

Catalog No	YP-mAb-03100
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CIB3
Immunogen	The antiserum was produced against synthesized peptide derived from human CIB3. AA range:138-187
Specificity	CIB3 Monoclonal Antibody detects endogenous levels of CIB3 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C7314M21Rik;Calcium and integrin binding family member 3;DNA dependent protein kinase catalytic subunit interacting protein 3;Gm117;Kinase interacting protein 3;KIP 3;KIP3;MGC13845;MGC142151;MGC96922;CIB3;Calcium and integrin-binding family member 3;Kinase-interacting protein 3
Observed Band	32kD
Calculated Molecular Weight	32kD
Cell Pathway	
Tissue Specificity	PCR rescued clones,
Function	similarity:Contains 3 EF-hand domains.,
Background	This gene product shares a high degree of sequence similarity with DNA-dependent protein kinase catalytic subunit-interacting protein 2 in human and mouse, and like them may bind the catalytic subunit of DNA-dependent protein kinases. The exact function of this gene is not known. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2014],


**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using CIB3 Monoclonal Antibody