



CLIP3 Monoclonal Antibody

Catalog No	YP-mAb-03103
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CLIP3
Immunogen	The antiserum was produced against synthesized peptide derived from human CLIP3. AA range:361-410
Specificity	CLIP3 Monoclonal Antibody detects endogenous levels of CLIP3 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	155P14Rik;AI844915;CAP GLY domain containing linker protein 3;CAP-Gly domain-containing linker protein 3;CLIP 17 related 59 kDa protein;CLIP-17-related 59 kDa protein;clip3;CLIP3_HUMAN;CLIPR 59;CLIPR-59;CLIPR59;Cytoplasmic linker protein 17 related 59 kDa protein;Cytoplasmic linker protein 17-related 59 kDa protein;Restin like 1;RSNL1;Cytoplasmic linker protein 170-related 59 kDa protein;CLIP-170-related 59 kDa protein
Observed Band	60kD
Calculated Molecular Weight	60kD
Cell Pathway	Cell membrane ; Lipid-anchor . Cytoplasm . Golgi apparatus, Golgi stack . Localized to Golgi stacks as well as on tubulovesicular elements juxtaposed to Golgi cisternae.
Tissue Specificity	Brain,Skin,Uterus,
Function	domain:Microtubule association is inhibited by the ANK repeats and the Golgi localization region (GoLD).,function:Functions as a cytoplasmic linker protein. Involved in TGN-endosome dynamics.,miscellaneous:The N-terminal half is dispensable for proper Golgi targeting, whereas the GoLD region is required.,similarity:Contains 2 CAP-Gly domains.,similarity:Contains 3 ANK repeats.,



subcellular location:Localized to Golgi stacks as well as on tubulovesicular elements juxtaposed to Golgi cisternae.,subunit:Homodimer.,

Background

This gene encodes a member of the cytoplasmic linker protein 170 family. Members of this protein family contain a cytoskeleton-associated protein glycine-rich domain and mediate the interaction of microtubules with cellular organelles. The encoded protein plays a role in T cell apoptosis by facilitating the association of tubulin and the lipid raft ganglioside GD3. The encoded protein also functions as a scaffold protein mediating membrane localization of phosphorylated protein kinase B. Alternatively spliced transcript variants have been observed for this gene. [provided by RefSeq, Dec 2010],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using CLIP3 Monoclonal Antibody