



GCP4 Monoclonal Antibody

Catalog No	YP-mAb-03137
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TUBGCP4
Immunogen	The antiserum was produced against synthesized peptide derived from human TUBGCP4. AA range:191-240
Specificity	GCP4 Monoclonal Antibody detects endogenous levels of GCP4 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	76P;GCP 4;GCP4;h76p;hGCP4;hGrip76;TUBGCP4;Gamma-tubulin complex component 4;GCP-4;Gamma-ring complex protein 76 kDa
Observed Band	76kD
Calculated Molecular Weight	76kD
Cell Pathway	Cytoplasm, cytoskeleton, microtubule organizing center, centrosome.
Tissue Specificity	Ubiquitously expressed.
Function	function:Gamma-tubulin complex is necessary for microtubule nucleation at the centrosome.,similarity:Belongs to the GCP family.,subunit:Gamma-tubulin complex is composed of gamma-tubulin, GCP2, GCP3, GCP4, GCP5 and GCP6.,tissue specificity:Ubiquitously expressed.,
Background	This gene encodes a component of the gamma-tubulin ring complex, which is required for microtubule nucleation. In mammalian cells, the protein localizes to centrosomes in association with gamma-tubulin. Crystal structure analysis revealed a structure composed of five helical bundles arranged around conserved hydrophobic cores. An exposed surface area located in the C-terminal domain is essential and sufficient for direct binding to gamma-tubulin.



Mutations in this gene that alter microtubule organization are associated with microcephaly and chorioretinopathy. Alternative splicing results in multiple transcript variants. [provided by RefSeq, May 2015],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using GCP4 Monoclonal Antibody