



MYH14 Monoclonal Antibody

Catalog No	YP-mAb-03155
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MYH14
Immunogen	The antiserum was produced against synthesized peptide derived from human MYH14. AA range:1051-1100
Specificity	MYH14 Monoclonal Antibody detects endogenous levels of MYH14 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Myosin 14;myosin;MHC16;KIAA2034;FP17425;MYH14;Myosin-14;Myosin heavy chain 14;Myosin heavy chain;non-muscle IIc;Non-muscle myosin heavy chain IIc;NMHC II-C
Observed Band	228kD
Calculated Molecular Weight	228kD
Cell Pathway	stress fiber,cytosol,brush border,membrane,myosin complex,myosin II complex, axon,growth cone,actomyosin,myelin sheath,extracellular exosome,myosin II filament,
Tissue Specificity	High levels of expression are found in brain (highest in corpus callosum), heart, kidney, liver, lung, small intestine, colon and skeletal muscle. Expression is low in organs composed mainly of smooth muscle, such as aorta, uterus and urinary bladder. No detectable expression is found in thymus, spleen, placenta and lymphocytes.
Function	disease:Defects in MYH14 are the cause of non-syndromic sensorineural deafness autosomal dominant type 4 (DFNA4) [MIM:600652]. DFNA4 is a form of sensorineural hearing loss. Sensorineural deafness results from damage to the neural receptors of the inner ear, the nerve pathways to the brain, or the area of the brain that receives sound information.,domain:The rodlike tail



sequence is highly repetitive, showing cycles of a 28-residue repeat pattern composed of 4 heptapeptides, characteristic for alpha-helical coiled coils., function:Cellular myosin that appears to play a role in cytokinesis, cell shape, and specialized functions such as secretion and capping.,sequence caution:Translation N-terminally extended.,similarity:Contains 1 IQ domain., similarity:Contains 1 myosin head-like domain.,subunit:Myosin is a hexameric protein that consists of 2 heavy chain subunits (MHC), 2 alkali light cha

Background

This gene encodes a member of the myosin superfamily. The protein represents a conventional non-muscle myosin; it should not be confused with the unconventional myosin-14 (MYO14). Myosins are actin-dependent motor proteins with diverse functions including regulation of cytokinesis, cell motility, and cell polarity. Mutations in this gene result in one form of autosomal dominant hearing impairment. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2011],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using MYH14 Monoclonal Antibody