



POTE14/22 Monoclonal Antibody

Catalog No	YP-mAb-03170
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	POTEG
Immunogen	The antiserum was produced against synthesized peptide derived from human POTEG. AA range:401-450
Specificity	POTE14/22 Monoclonal Antibody detects endogenous levels of POTE14/22 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ANKRD26-like family C member 2;ovary;POTE ankyrin domain family member G;POTE-14;POTEG;POTEG_HUMAN;Prostate;testis-expressed protein on chromosome 14;ACT;ANKRD26 like family C member ENSP34942;P74P;POTE ankyrin domain family member M;POTE14beta;POTEM;POTEM_HUMAN;Prostate specific P74P;Putative POTE ankyrin domain family member M;POTE14/22;A26C2;POTE14;POTEH;A26C3;ACTBL1;POTE22;POTE ankyrin domain family member H;ANKRD26-like; testis-expressed protein on chromosome 14
Observed Band	67kD
Calculated Molecular Weight	67kD
Cell Pathway	
Tissue Specificity	Prostate,
Function	similarity:Belongs to the POTE family.,similarity:Contains 5 ANK repeats.,
Background	similarity:Belongs to the POTE family.,similarity:Contains 5 ANK repeats.,




matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using POTE14/22 Monoclonal Antibody