



POTE-14/22 Monoclonal Antibody

Catalog No	YP-mAb-03171
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	POTEG/POTEH
Immunogen	The antiserum was produced against synthesized peptide derived from human A26C2/3. AA range:178-227
Specificity	POTE-14/22 Monoclonal Antibody detects endogenous levels of POTE-14/22 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ANKRD26-like family C member 2;ovary;POTE ankyrin domain family member G;POTE-14;POTEG;POTEG_HUMAN;Prostate;testis-expressed protein on chromosome 14;POTE-14/22;POTEG/POTEH;POTE ankyrin domain family member G/H;POTEH;A26C3;ACTBL1;POTE22;POTE ankyrin domain family member H;ANKRD26-like family C member 3;testis-expressed protein on chromosome 22;POTE-22;A26C2;POTE14;ANKRD26-like; testis-expressed protein on chromosome 22
Observed Band	61kD
Calculated Molecular Weight	61kD
Cell Pathway	integral component of membrane,
Tissue Specificity	Prostate,
Function	similarity:Belongs to the POTE family.,similarity:Contains 7 ANK repeats.,
Background	similarity:Belongs to the POTE family.,similarity:Contains 7 ANK repeats.,


**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using POTE-14/22 Monoclonal Antibody