




POTE15 Monoclonal Antibody

Catalog No	YP-mAb-03172
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	POTEB
Immunogen	The antiserum was produced against synthesized peptide derived from human A26B1. AA range:481-530
Specificity	POTE15 Monoclonal Antibody detects endogenous levels of POTE15 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	POTEB;A26B1;ANKRD26 like family B member 1;POTE 15;POTE ankyrin domain family member B3;POTE15;POTE ankyrin domain family member B; ANKRD26-like family B member 1;Prostate;ovary;testis-expressed protein on chromosome 15;POTE-15; testis-expressed protein on chromosome 15
Observed Band	60kD
Calculated Molecular Weight	60kD
Cell Pathway	
Tissue Specificity	
Function	
Background	POTEB (POTE Ankyrin Domain Family Member B) is a Protein Coding gene. An important paralog of this gene is ANKRD18A.
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images**Mouse-kidney 3T3**

Western Blot analysis of various cells using POTE15 Monoclonal Antibody