




S-100A7L2 Monoclonal Antibody

Catalog No	YP-mAb-03182
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	S100A7L2
Immunogen	The antiserum was produced against synthesized peptide derived from human S100A7L2. AA range:38-87
Specificity	S-100A7L2 Monoclonal Antibody detects endogenous levels of S-100A7L2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Protein S1 A7 like 2;S1 calcium binding protein A7 like 2;S1A7B;S1A7L2;S-100A7L2;S100A7L2;Protein S100-A7-like 2;S100A7B;S100 calcium-binding protein A7-like 2
Observed Band	11kD
Calculated Molecular Weight	11kD
Cell Pathway	
Tissue Specificity	
Function	similarity:Belongs to the S-100 family.,similarity:Contains 2 EF-hand domains.,
Background	similarity:Belongs to the S-100 family.,similarity:Contains 2 EF-hand domains.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using S-100A7L2
Monoclonal Antibody