



S-100Z Monoclonal Antibody

Catalog No	YP-mAb-03184
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	S100Z
Immunogen	The antiserum was produced against synthesized peptide derived from human S100Z. AA range:2-51
Specificity	S-100Z Monoclonal Antibody detects endogenous levels of S-100Z protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Gm625;Protein S1 Z;Protein S1-Z;S1 calcium binding protein Z;S1 calcium binding protein zeta;S1 calcium-binding protein Z;S1 zeta;S1Z;S1Z_HUMAN;S-100Z;S100Z;Protein S100-Z;S100 calcium-binding protein Z
Observed Band	12kD
Calculated Molecular Weight	12kD
Cell Pathway	
Tissue Specificity	Highest level of expression in spleen and leukocytes.
Function	miscellaneous:This protein binds two calcium ions.,similarity:Belongs to the S-100 family.,similarity:Contains 2 EF-hand domains.,subunit:Homodimer. Interacts with S100P.,tissue specificity:Highest level of expression in spleen and leukocytes.,
Background	Members of the S100 protein family contain 2 calcium-binding EF-hands and exhibit cell-type specific expression patterns. For additional background information on S100 proteins, see MIM 114085.[supplied by OMIM, Mar 2008],


**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using S-100Z Monoclonal Antibody