



WAVE3 Monoclonal Antibody

Catalog No	YP-mAb-03216
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	WASF3
Immunogen	The antiserum was produced against synthesized peptide derived from human WASF3. AA range:151-200
Specificity	WAVE3 Monoclonal Antibody detects endogenous levels of WAVE3 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Brush 1;KIAA0900;Protein WAVE 3;SCAR3;WAS protein family;member 3; WASF3;WASP family protein member 3;WAVE3;Wiskott-Aldrich syndrome protein family member 3;Protein WAVE-3;Verprolin homology domain-containing protein 3
Observed Band	55kD
Calculated Molecular Weight	55kD
Cell Pathway	Cytoplasm, cytoskeleton.
Tissue Specificity	Expressed in ovary and brain.
Function	domain: Binds the Arp2/3 complex through the C-terminal region and actin through verprolin homology (VPH) domain.,function: Downstream effector molecules involved in the transmission of signals from tyrosine kinase receptors and small GTPases to the actin cytoskeleton.,similarity: Belongs to the SCAR/ WAVE family.,similarity: Contains 1 WH2 domain.,subunit: Binds actin and the Arp2/3 complex.,tissue specificity: Expressed in ovary and brain.,



Background

This gene encodes a member of the Wiskott-Aldrich syndrome protein family. The gene product is a protein that forms a multiprotein complex that links receptor kinases and actin. Binding to actin occurs through a C-terminal verprolin homology domain in all family members. The multiprotein complex serves to transduce signals that involve changes in cell shape, motility or function. A pseudogene of this gene have been defined on chromosome 6. Alternative splicing results in multiple transcript variants [provided by RefSeq, May 2014],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using WAVE3 Monoclonal Antibody