



ARP-1 Monoclonal Antibody

Catalog No	YP-mAb-03307
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	NR2F2
Immunogen	The antiserum was produced against synthesized peptide derived from human COT2. AA range:1-50
Specificity	ARP-1 Monoclonal Antibody detects endogenous levels of ARP-1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Apolipoprotein A-I regulatory protein 1;Apolipoprotein AI regulatory protein 1; ARP-1;ARP1;COT2_HUMAN;COUP TF2;COUP transcription factor 2;COUP transcription factor II;COUP-TF II;COUP-TF2;NR2F2;Nuclear receptor subfamily 2 group F member 2;TFCOUP2
Observed Band	45kD
Calculated Molecular Weight	45kD
Cell Pathway	Nucleus.
Tissue Specificity	Ubiquitous. Expressed in the stromal cells of developing fetal ovaries (PubMed:29478779).
Function	function:Regulation of the apolipoprotein A-I gene transcription. Binds to DNA site A.,similarity:Belongs to the nuclear hormone receptor family., similarity:Belongs to the nuclear hormone receptor family. NR2 subfamily., similarity:Contains 1 nuclear receptor DNA-binding domain.,subunit:Interacts with SQSTM1 (By similarity). Binds DNA as a homodimer.,tissue specificity:Found in all tissues.,



Background

This gene encodes a member of the steroid thyroid hormone superfamily of nuclear receptors. The encoded protein is a ligand inducible transcription factor that is involved in the regulation of many different genes. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Mar 2010],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using ARP-1 Monoclonal Antibody