



PGRMC2 Monoclonal Antibody

Catalog No	YP-mAb-03323
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PGRMC2
Immunogen	The antiserum was produced against synthesized peptide derived from human PGRMC2. AA range:174-223
Specificity	PGRMC2 Monoclonal Antibody detects endogenous levels of PGRMC2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DG6;PMBP;PGRMC2;Membrane-associated progesterone receptor component 2;Progesterone membrane-binding protein;Steroid receptor protein DG6
Observed Band	25kD
Calculated Molecular Weight	25kD
Cell Pathway	Membrane ; Single-pass membrane protein . Nucleus envelope . Endoplasmic reticulum .
Tissue Specificity	Expressed by endometrial glands and stroma (at protein level).
Function	domain:The cytochrome b5 heme-binding domain lacks the conserved iron-binding His residues at positions 137 and 161.,function:Receptor for steroids ., similarity:Belongs to the cytochrome b5 family. MAPR subfamily., similarity:Contains 1 cytochrome b5 heme-binding domain.,
Background	domain:The cytochrome b5 heme-binding domain lacks the conserved iron-binding His residues at positions 137 and 161.,function:Receptor for steroids ., similarity:Belongs to the cytochrome b5 family. MAPR subfamily., similarity:Contains 1 cytochrome b5 heme-binding domain.,




matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using PGRMC2 Monoclonal Antibody