



4.1G Monoclonal Antibody

Catalog No	YP-mAb-03665
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	EPB41L2
Immunogen	Synthesized peptide derived from 4.1G . at AA range: 530-610
Specificity	4.1G Monoclonal Antibody detects endogenous levels of 4.1G protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	4.1 G;4.1G;AW555191;Band 4.1 like protein 2;Band 4.1-like protein 2; D1ErtD398e;DKFZp781D1972;DKFZp781H1755;E41L2_HUMAN;Epb41l2; Erythrocyte membrane protein band 4.1 like 2;Erythrocyte membrane protein band 4.1 like protein 2;Generally expressed protein 4.1;NBL2;OTTHUMP17197; OTTHUMP4262
Observed Band	112kD
Calculated Molecular Weight	112kD
Cell Pathway	Cytoplasm, cytoskeleton . Cytoplasm, cell cortex . Cell membrane .
Tissue Specificity	Widely expressed.
Function	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR., similarity:Contains 1 FERM domain.,subunit:The CTD domain interacts with FKBP2 (By similarity). Interacts with FCGR1A.,tissue specificity:Widely expressed.,
Background	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR., similarity:Contains 1 FERM domain.,subunit:The CTD domain interacts with FKBP2 (By similarity). Interacts with FCGR1A.,tissue specificity:Widely expressed.,




matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using 4.1G Monoclonal Antibody