



AKAP 14 Monoclonal Antibody

Catalog No	YP-mAb-03689
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	AKAP14
Immunogen	The antiserum was produced against synthesized peptide derived from human AKAP14. AA range:1-50
Specificity	AKAP 14 Monoclonal Antibody detects endogenous levels of AKAP 14 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	A kinase anchor protein 14;A kinase anchor protein 28 kDa;AKAP 14;AKAP 28; AKAP14;AKAP28;A-kinase anchor protein 14;AKAP-14;A-kinase anchor protein 28 kDa;Protein kinase A-anchoring protein 14;PRKA14
Observed Band	23kD
Calculated Molecular Weight	23kD
Cell Pathway	Cytoplasm .
Tissue Specificity	Present in cilia (at protein level). Expressed in tissues containing axoneme-based organelles (cilia and/or flagella): trachea and testis. Highly expressed in airway cilia.
Function	caution:The sequence shown here is derived from an Ensembl automatic analysis pipeline and should be considered as preliminary data.,function:Binds to type II regulatory subunits of protein kinase A and anchors/targets them., subunit:Binds to type II regulatory subunits (RII).,tissue specificity:Present in cilia (at protein level). Expressed in tissues containing axoneme-based organelles (cilia and/or flagella): trachea and testis. Highly expressed in airway cilia.,



Background

The A-kinase anchor proteins (AKAPs) are a group of structurally diverse proteins, which have the common function of binding to the regulatory subunit of protein kinase A (PKA) and confining the holoenzyme to discrete locations within the cell. This gene encodes a member of the AKAP family. The protein anchors PKA in ciliary axonemes and, in this way, may play a role in regulating ciliary beat frequency. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using AKAP 14 Monoclonal Antibody