




Annexin X Monoclonal Antibody

Catalog No	YP-mAb-03705
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ANXA10
Immunogen	The antiserum was produced against synthesized peptide derived from human ANXA10. AA range:261-310
Specificity	Annexin X Monoclonal Antibody detects endogenous levels of Annexin X protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Annexin A10;Annexin 10;Annexin 14;ANX10;ANX14;Annexin X;ANXA10; Annexin-10;Annexin-14
Observed Band	38kD
Calculated Molecular Weight	38kD
Cell Pathway	mitochondrion,
Tissue Specificity	Colon,Lung,
Function	similarity:Belongs to the annexin family.,similarity:Contains 4 annexin repeats.,
Background	This gene encodes a member of the annexin family. Members of this calcium-dependent phospholipid-binding protein family play a role in the regulation of cellular growth and in signal transduction pathways. The function of this gene has not yet been determined. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using Annexin X Monoclonal Antibody