




ARMC6 Monoclonal Antibody

Catalog No	YP-mAb-03713
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ARMC6
Immunogen	The antiserum was produced against synthesized peptide derived from human ARMC6. AA range:431-480
Specificity	ARMC6 Monoclonal Antibody detects endogenous levels of ARMC6 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Armadillo repeat-containing protein 6;R30923_1;ARMC6
Observed Band	54kD
Calculated Molecular Weight	54kD
Cell Pathway	
Tissue Specificity	Cervix,Hippocampus,Melanoma,Skin,
Function	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR., similarity:Belongs to the ARMC6 family.,similarity:Contains 4 ARM repeats.,
Background	The function of this gene's protein product has not been determined. A related protein in mouse suggests that this protein has a conserved function. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Nov 2010],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using ARMC6 Monoclonal Antibody