




Arrdc1 Monoclonal Antibody

Catalog No	YP-mAb-03718
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ARRDC1
Immunogen	The antiserum was produced against synthesized peptide derived from human ARRD1. AA range:167-216
Specificity	Arrdc1 Monoclonal Antibody detects endogenous levels of Arrdc1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Alpha-arrestin 1;arrestin domain containing 1;Arrestin domain-containing protein 1;Arrdc1
Observed Band	48kD
Calculated Molecular Weight	48kD
Cell Pathway	Cell membrane . Also found in extracellular vesicles different from exosomes. .
Tissue Specificity	Blood,Brain,
Function	similarity:Belongs to the arrestin family.,
Background	similarity:Belongs to the arrestin family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using Arrdc1
Monoclonal Antibody