




ARRDC2 Monoclonal Antibody

Catalog No	YP-mAb-03719
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ARRDC2
Immunogen	The antiserum was produced against synthesized peptide derived from human ARRDC2. AA range:141-190
Specificity	ARRDC2 Monoclonal Antibody detects endogenous levels of ARRDC2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ARRDC2;Arrestin domain-containing protein 2;PP273;PP2703
Observed Band	44kD
Calculated Molecular Weight	44kD
Cell Pathway	plasma membrane,cytoplasmic vesicle,
Tissue Specificity	Brain,Synovial cell,
Function	similarity:Belongs to the arrestin family.,
Background	similarity:Belongs to the arrestin family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using ARRDC2 Monoclonal Antibody