




C1qL2 Monoclonal Antibody

Catalog No	YP-mAb-03742
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	C1QL2
Immunogen	The antiserum was produced against synthesized peptide derived from human C1QL2. AA range:231-280
Specificity	C1qL2 Monoclonal Antibody detects endogenous levels of C1qL2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C1q and tumor necrosis factor-related protein 10;C1q/TNF-related protein 10; Complement C1q-like protein 2;CTRP10;C1qL2
Observed Band	32-36kD
Calculated Molecular Weight	32-36kD
Cell Pathway	Secreted .
Tissue Specificity	
Function	similarity:Contains 1 C1q domain.,similarity:Contains 1 collagen-like domain.,
Background	similarity:Contains 1 C1q domain.,similarity:Contains 1 collagen-like domain.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using C1qL2
Monoclonal Antibody