



CD172b Monoclonal Antibody

Catalog No	YP-mAb-03749
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SIRPB1
Immunogen	The antiserum was produced against synthesized peptide derived from human SIRPB1. AA range:231-280
Specificity	CD172b Monoclonal Antibody detects endogenous levels of CD172b protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SIRPB1;SIRP Beta 1;SIRP-beta-1;Signal-regulatory protein beta-1;CD172b; CD172 antigen-like family member B;CD antigen CD172b
Observed Band	48kD
Calculated Molecular Weight	48kD
Cell Pathway	Cell membrane ; Single-pass type I membrane protein .
Tissue Specificity	Detected in monocytes and dendritic cells.
Function	function:Immunoglobulin-like cell surface receptor involved in the negative regulation of receptor tyrosine kinase-coupled signaling processes., function:Immunoglobulin-like cell surface receptor involved in the negative regulation of receptor tyrosine kinase-coupled signaling processes. Participates also in the recruitment of tyrosine kinase SYK.,PTM:N-glycosylated., similarity:Contains 1 Ig-like V-type (immunoglobulin-like) domain., similarity:Contains 2 Ig-like C1-type (immunoglobulin-like) domains., subunit:Interacts with TYROBP. This interaction results in the recruitment of SYK.,tissue specificity:Detected in monocytes and dendritic cells.,



Background

The protein encoded by this gene is a member of the signal-regulatory-protein (SIRP) family, and also belongs to the immunoglobulin superfamily. SIRP family members are receptor-type transmembrane glycoproteins known to be involved in the negative regulation of receptor tyrosine kinase-coupled signaling processes. This protein was found to interact with TYROBP/DAP12, a protein bearing immunoreceptor tyrosine-based activation motifs. This protein was also reported to participate in the recruitment of tyrosine kinase SYK. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Feb 2009],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using CD172b Monoclonal Antibody