



CEP152 Monoclonal Antibody

Catalog No	YP-mAb-03757
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CEP152
Immunogen	The antiserum was produced against synthesized peptide derived from human CEP152. AA range:841-890
Specificity	CEP152 Monoclonal Antibody detects endogenous levels of CEP152 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KIAA0912;centrosomal protein 152kDa;CEP 152;CEP152;Centrosomal protein of 152 kDa
Observed Band	140kD
Calculated Molecular Weight	140kD
Cell Pathway	Cytoplasm, cytoskeleton, microtubule organizing center, centrosome . Cytoplasm, cytoskeleton, microtubule organizing center, centrosome, centriole . Colocalizes with CDK5RAP2, WDR62 and CEP63 in a discrete ring around the proximal end of the parental centriole. At this site, a cohesive structure is predicted to engage parental centrioles and procentrioles (PubMed:21983783, PubMed:26297806). Localizes to the deuterosome (By similarity). Localizes to pericentriolar material (PCM) (PubMed:26337392). .
Tissue Specificity	Bladder,Brain,Heart,Lung,Testis,
Function	sequence caution:Contaminating sequence. Potential poly-A sequence.,
Background	centrosomal protein 152(CEP152) Homo sapiens This gene encodes a protein that is thought to be involved with centrosome function. Mutations in this gene have been associated with primary microcephaly (MCPH4). Alternative splicing



results in multiple transcript variants. [provided by RefSeq, Aug 2010],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using CEP152 Monoclonal Antibody