




CEP70 Monoclonal Antibody

Catalog No	YP-mAb-03764
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CEP70
Immunogen	The antiserum was produced against synthesized peptide derived from human CEP70. AA range:241-290
Specificity	CEP70 Monoclonal Antibody detects endogenous levels of CEP70 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	p10-binding protein;BITE;CEP70;Centrosomal protein of 70 kDa
Observed Band	70kD
Calculated Molecular Weight	70kD
Cell Pathway	Cytoplasm, cytoskeleton, microtubule organizing center, centrosome . Localized at the center of the radial microtubule array in interphase and at the spindle poles during various stages of mitosis.
Tissue Specificity	Chronic myeloid leukemia cell,Skin,Teratocarcinoma,Testis,
Function	similarity:Contains 1 TPR repeat.,
Background	similarity:Contains 1 TPR repeat.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using CEP70 Monoclonal Antibody