



CEP78 Monoclonal Antibody

Catalog No	YP-mAb-03766
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CEP78
Immunogen	The antiserum was produced against synthesized peptide derived from human CEP78. AA range:301-350
Specificity	CEP78 Monoclonal Antibody detects endogenous levels of CEP78 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C9orf81;centrosomal protein 78kDa;Centrosomal protein of 78 kDa;CEP78;IP63
Observed Band	80kD
Calculated Molecular Weight	80kD
Cell Pathway	Cytoplasm, cytoskeleton, microtubule organizing center, centrosome . Cytoplasm, cytoskeleton, microtubule organizing center, centrosome, centriole . Cytoplasm, cytoskeleton, cilium basal body . Mainly localizes at the centriolar wall, but also found in the pericentriolar material (PubMed:27246242). Expressed in photoreceptor inner segment (PubMed:27588452). .
Tissue Specificity	Widely expressed (PubMed:27588451, PubMed:27588452). Expressed in different retinal cell types with higher expression in cone compared to rod cells (at protein level) (PubMed:27588452).
Function	similarity:Belongs to the CEP78 family.,similarity:Contains 3 LRR (leucine-rich) repeats.,
Background	similarity:Belongs to the CEP78 family.,similarity:Contains 3 LRR (leucine-rich) repeats.,




matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using CEP78 Monoclonal Antibody