



CHST2 Monoclonal Antibody

Catalog No	YP-mAb-03772
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CHST2
Immunogen	The antiserum was produced against synthesized peptide derived from human CHST2. AA range:1-50
Specificity	CHST2 Monoclonal Antibody detects endogenous levels of CHST2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C6ST;EC:2.8.2.-;Galactose/N-acetylglucosamine/N-acetylglucosamine 6-O-sulfotransferase 2;GlcNAc6ST 1;GlcNAc6ST-1;CHST2;Carbohydrate sulfotransferase 2;GN6ST;GST-2;N-acetylglucosamine 6-O-sulfotransferase 1; Gn6ST-1
Observed Band	58kD
Calculated Molecular Weight	58kD
Cell Pathway	Golgi apparatus, trans-Golgi network membrane ; Single-pass type II membrane protein .
Tissue Specificity	Widely expressed. Highly expressed in bone marrow, peripheral blood leukocytes, spleen, brain, spinal cord, ovary and placenta. Expressed by high endothelial cells (HEVs) and leukocytes.
Function	caution:It is uncertain whether Met-1 or Met-48 is the initiator.,function:Catalyzes the transfer of sulfate to position 6 of non-reducing N-acetylglucosamine (GlcNAc) residues within keratan-like structures on N-linked glycans and within mucin-associated glycans that can ultimately serve as L-selectin ligands. L-selectin ligands are present in high endothelial cells (HEVs) and play a central role in lymphocyte homing at sites of inflammation. Participates in biosynthesis of L-selectin ligand sialyl 6-sulfo Lewis X and in lymphocyte homing to Peyer



patches. Has no activity toward O-linked sugars. Its substrate specificity may be influenced by its subcellular location. Sulfates GlcNAc residues at terminal, non-reducing ends of oligosaccharide chains.,induction:Up-regulated upon cytokine activation.,online information:GlycoGene database,similarity:Belongs to the sulfotransferase 1 family.

Background

This locus encodes a sulfotransferase protein. The encoded enzyme catalyzes the sulfation of a nonreducing N-acetylglucosamine residue, and may play a role in biosynthesis of 6-sulfosialyl Lewis X antigen. [provided by RefSeq, Aug 2011],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using CHST2 Monoclonal Antibody